

Introduction:

High-grade serous carcinoma (HGSC) is the most common subtype of ovarian cancer, with a high mortality rate, even after successful first-line treatment (debulking surgery and chemotherapy). Although therapeutic options for targeted therapy are rapidly expanding, identification of patients who respond to these therapies remains a challenge. In this study, we used a new technology that measures functional signal transduction pathway (STP) activity to identify HGSCs with a potential clinical target by comparing STP activity in HGSC and in normal Fallopian tube epithelium (FTE); the tissue of origin of most HGSCs.

Methods:

Samples:

50 primary tissue samples from postmenopausal patients with advanced stage HGSC, PRIOR to start of chemotherapy and 9 postmenopausal morphologically normal FTE samples (by laser capture microdissection).

Pathway analysis:

STP activity was measured using the G7 (qPCR) OncoSIGNal test (InnoSIGN BV, The Netherlands). Measured pathways: androgen receptor (AR), estrogen receptor (ER), PI3K, MAPK, Hedgehog (HH), TGFβ and Notch.

Statistics:


Differences in STP activity between groups were compared with Mann-Whitney U tests. Cut-off value for aberrant STP activity was defined as the 95th percentile of STP activity measured in FTE.

Results:

HGSC group (n=50) as compared to FTE (n=9):

Lower median ER ($p < 0.001$) pathway activity.

Higher median PI3K ($p < 0.01$), HH ($p < 0.001$) and TGFβ pathway activity ($p = 0.01$).



Results:

Individual HGSC samples:

Aberrant increased activity was identified for the MAPK (n = 11), PI3K (n = 22), HH (n = 27) and TGFβ (n = 19) pathways.

Frequently observed combinations of aberrant STP activity were HH/TGFβ (n = 15) and HH/PI3K (n = 9).

ID	OSST19791	OSST19809	OSST19810	OSST19794	OSST18609	OSST18664	OSST19779	OSST19787	OSST19817	OSST19808	OSST19796	OSST19786	OSST19801	OSST19819	OSST19806	OSST19781	OSST19793	OSST19814	OSST19821	OSST19788	OSST19792	OSST18614	OSST18612	OSST19798	OSST19783	OSST19795	OSST19800	OSST19804	OSST19811	OSST19818	OSST19789	OSST19820	OSST19822	OSST19816	OSST19785	OSST19805	OSST19824	OSST18663	OSST19797	OSST19802	OSST19803	OSST19812	OSST19813	OSST19823	OSST19784	OSST19780	OSST19782	OSST19790	OSST19815							
HH	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓		
PI3K																																																								
TGFβ																																																								
MAPK																																																								

Aberrant increased STP activity indicated per HGSC patient by color-coding. Distinct aberrant STP profiles among the HGSC cohort.

Conclusion

- In 92% of analyzed HGSC samples a clinical target could be identified using STP analysis.
- Differentiation between normal and aberrant STP activity could have clinical utility in the selection of HGSC patients for targeted therapy.

Next step:

A prospective study in ovarian cancer patients (STAPOVER) has been designed to demonstrate clinical utility (NCT03458221). Additional FTE samples will be measured to increase sample size.

