


oncosignal

Laboratory Developed Test

Instructions for submission



InnoSIGN's LDT for patients with invasive breast cancer

Instructions for submission of pathology material to InnoSIGN for OncoSIGNal testing / evaluation of invasive breast cancers.

1. The test requisition form will be provided to the Pathology Department directly from the Medical Oncologist.
2. Please obtain a complete Pathology report on the selected patient which represents the positive biopsy or positive surgical excision specimen.
3. Please show the Pathology report to the Pathologist who signed out the case and ask him/her to review the Pathology report and select the specimen part in the report that contains the most tumor.
4. Please pull the slides that are associated with that specimen part and ask the Pathologist to pick the **one slide that is the best representation of this patient's invasive breast cancer.**
5. Please locate the Formalin Fixed Paraffin Embedded tissue block associated with this histology slide. Mail the selected histology slide and add either **the associated tissue block** or 5 slides of 10µm prepared from associated tissue block, with at least 10 mm² of tumor tissue per slide, together with **the completed test requisition form** and **the pathology report** to InnoSIGN Laboratory in Mason, Ohio.
6. Please use the FedEx shipping label provided by InnoSIGN.
7. Address:
InnoSIGN
5155 Financial Way, Suite 17
Mason, OH, 45040
8. Phone number: 800-686-0844 (toll-free)
9. E-mail: customersupport@innosignbio.com.

The histology slide and associated tissue block will be returned to your department once the testing is completed.

We appreciate the effort of all involved with this process on behalf of this patient.

Without the help of the Pathologist and the support of all the staff in the Department of Pathology - no such testing could take place. We are grateful for your help.